

# Advanced Guide to Ventilators and Ventilator Settings

## Overview of Mechanical Ventilation

Mechanical ventilation is a cornerstone of critical care, providing life-sustaining support for patients with acute respiratory failure (ARF), acute respiratory distress syndrome (ARDS), or other conditions compromising gas exchange or ventilatory mechanics. It delivers positive pressure ventilation via an endotracheal tube or tracheostomy, optimizing oxygenation and CO<sub>2</sub> elimination while minimizing ventilator-induced lung injury (VILI). This advanced guide is intended for clinicians with a strong medical background, offering a detailed exploration of ventilator physiology, settings, troubleshooting, clinical applications, and optimization strategies, with tables and clinical scenarios for practical application.

## Advanced Physiology of Mechanical Ventilation

### Physiological Principles:

- Gas Exchange:
  - **Oxygenation:** Determined by FiO<sub>2</sub>, PEEP, and mean airway pressure (MAP).  $MAP = (P_{peak} \times T_i + PEEP \times T_e) / (T_i + T_e)$ , where  $T_i$  is inspiratory time and  $T_e$  is expiratory time. Higher MAP improves PaO<sub>2</sub> but risks barotrauma.
  - **Ventilation:** Governed by minute ventilation ( $VE = VT \times RR$ ). Alveolar ventilation ( $VA = [VT - V_d] \times RR$ , where  $V_d$  is dead space) directly influences PaCO<sub>2</sub>. Dead space fraction ( $V_d/VT$ ) increases in ARDS, PE, or hypovolemia, reducing effective ventilation.
- Lung Mechanics
  - **Compliance (C):**  $\Delta V / \Delta P$  (normal ~100 mL/cmH<sub>2</sub>O); decreased in ARDS, pulmonary edema, or fibrosis ( $C < 50$  mL/cmH<sub>2</sub>O). Static compliance ( $C_{st}$ ) is measured as  $VT / (P_{plat} - PEEP)$ .
  - **Resistance (R):**  $\Delta P / \Delta Flow$  (normal ~2-3 cmH<sub>2</sub>O/L/s); increased in bronchospasm, secretions, or ETT obstruction.  $R = (P_{peak} - P_{plat}) / Flow$ .
  - **Peak Pressure (P<sub>peak</sub>):** Reflects airway resistance and compliance;  $P_{peak} > 40$  cmH<sub>2</sub>O indicates obstruction or low compliance.
  - **Plateau Pressure (P<sub>plat</sub>):** Reflects alveolar pressure;  $P_{plat} > 30$  cmH<sub>2</sub>O risks barotrauma. Measured during an end-inspiratory pause (0.5-1 s hold).
- Ventilator-Induced Lung Injury (VILI):
  - **Barotrauma:** High  $P_{plat} (> 30 \text{ cmH}_2\text{O})$  causes alveolar overdistension.

- **Volutrauma:** High VT (>8 mL/kg IBW) overdistends alveoli, leading to inflammation.
- **Atelectrauma:** Low PEEP causes cyclic alveolar collapse/reopening, inducing shear stress.
- **Biotrauma:** Inflammatory cytokine release due to mechanical stress (e.g., TNF- $\alpha$ , IL-6).
  - **Lung-protective strategy:** VT 4-6 mL/kg IBW, Pplat <30 cmH<sub>2</sub>O, PEEP titrated to optimize oxygenation while minimizing overdistension.
- Oxygen Toxicity: Prolonged FiO<sub>2</sub> >60% risks oxidative stress, leading to absorption atelectasis, hyperoxic lung injury, and reactive oxygen species (ROS) production. Target SpO<sub>2</sub> 88-95% to minimize toxicity.
- Hemodynamic Effects:
  - Positive pressure increases intrathoracic pressure, reducing venous return and cardiac output (especially with high PEEP or MAP). RV afterload increases with high PEEP, potentially causing RV dysfunction (monitor with ECHO: RV dilation, septal flattening).

## Advanced Ventilator Settings: Indications, Optimization, and Rationale

### Modes of Ventilation:

- **Assist-Control (AC) Volume-Control (VC):**
  - Settings:
    - **VT:** 4-6 mL/kg IBW (ARDS), 6-8 mL/kg IBW (normal lungs).
    - **RR:** 12-20 breaths/min, titrate to PaCO<sub>2</sub> 35-45 mmHg (or permissive hypercapnia in ARDS, pH >7.20).
    - **FiO<sub>2</sub>:** Start at 100%, titrate to SpO<sub>2</sub> 88-95% (PaO<sub>2</sub> 60-100 mmHg).
    - **PEEP:** 5-15 cmH<sub>2</sub>O (higher in ARDS, titrate using PEEP-FiO<sub>2</sub> table per ARDSNet).
  - I:E ratio: 1:2 (normal), 1:1 (ARDS for oxygenation), 1:4 (obstructive diseases).
  - Indications: Severe ARF (e.g., ARDS, pneumonia), high WOB, or neuromuscular failure (e.g., myasthenia gravis).
  - Rationale: Provides full ventilatory support, ensures consistent VT delivery, reduces WOB. Patient-triggered breaths receive the set VT, preventing hypoventilation.

- **Assist-Control Pressure-Control (PC):**

- Settings:
  - **Inspiratory Pressure (P<sub>insp</sub>):** 15-25 cmH<sub>2</sub>O, titrate to VT 4-6 mL/kg IBW.
  - **RR:** 12-20 breaths/min.
  - **FiO<sub>2</sub>:** Titrate to SpO<sub>2</sub> 88-95%.
  - **PEEP:** 5-15 cmH<sub>2</sub>O.
  - **Ti:** 0.8-1.2 s (adjust for I:E ratio).
- Indications: ARDS, when precise pressure control is needed to limit P<sub>plat</sub> (<30 cmH<sub>2</sub>O).
- Rationale: Limits peak pressures, reducing barotrauma risk. VT varies with lung compliance, requiring frequent monitoring.

- **Synchronized Intermittent Mandatory Ventilation (SIMV):**

- Settings:
  - **VT:** 6-8 mL/kg IBW.
  - **RR:** 10-14 breaths/min (partial support).
  - **FiO<sub>2</sub>:** Titrate to SpO<sub>2</sub> 88-95%.
  - **PEEP:** 5-10 cmH<sub>2</sub>O.
  - **Pressure Support (PS):** 5-15 cmH<sub>2</sub>O for spontaneous breaths (target VT 4- 6 mL/kg).
- Indications: Weaning, patients with stable respiratory drive (e.g., post-surgery, resolving ARF).
- Rationale: Combines mandatory breaths with spontaneous breaths, reducing WOB while promoting patient effort. PS augments spontaneous breaths to prevent fatigue.

- **Pressure Support Ventilation (PSV):**

- Settings:
  - **PS:** 5-15 cmH<sub>2</sub>O (titrate to VT 4-6 mL/kg).
  - **FiO<sub>2</sub>:** 30-40% (weaning).
  - **PEEP:** 5 cmH<sub>2</sub>O.
- Indications: Spontaneous breathing trials (SBTs), weaning readiness assessment.
- Rationale: Supports spontaneous breaths, assesses patient ability to sustain ventilation without full support. Rapid shallow breathing index (RSBI = RR/VT in L) <105 predicts successful extubation.

- **Airway Pressure Release Ventilation (APRV):**

- Settings:
  - **P-high:** 20-30 cmH<sub>2</sub>O (titrate to oxygenation, MAP 15-20 cmH<sub>2</sub>O).
  - **P-low:** 0-5 cmH<sub>2</sub>O (optimize CO<sub>2</sub> clearance).
  - **T-high:** 4-6 s (prolonged high-pressure phase).

- **T-low:** 0.5-1 s (set to 50-75% of peak expiratory flow rate for optimal release).
  - **FiO2:** Titrate to SpO2 88-95%.
- Indications: Refractory hypoxemia in ARDS, severe atelectasis, or lung recruitment needs.
- Rationale: Maximizes alveolar recruitment (P-high, T-high), minimizes VILI by limiting cyclic collapse, allows spontaneous breathing to reduce sedation needs.
- **High-Frequency Oscillatory Ventilation (HFOV):**
  - Settings:
    - **Frequency:** 3-10 Hz (lower for larger patients).
    - **Mean Airway Pressure (mPaw):** 5-10 cmH2O above MAP on conventional ventilation.
    - **Amplitude (ΔP):** 60-90 cmH2O (titrate to chest wiggle to diaphragm).
    - **FiO2:** Titrate to SpO2 88-95%.
  - Indications: Rescue therapy for severe ARDS (PaO2/FiO2 <100 despite optimal conventional ventilation).
  - Rationale: Delivers very small VT (1-3 mL/kg) at high frequency, minimizing volutrauma while maintaining oxygenation via gas diffusion.

## Advanced Settings and Optimization:

---

- Tidal Volume (VT):
  - **ARDSNet Protocol:** VT 4-6 mL/kg IBW, Pplat <30 cmH2O to prevent volutrauma.
  - **Calculate IBW:** Men:  $50 + 2.3 \times [\text{height (in)} - 60]$ ; Women:  $45.5 + 2.3 \times [\text{height (in)} - 60]$ .
  - **Rationale:** Low VT reduces alveolar overdistension, decreases biotrauma (inflammatory cytokine release).
- Respiratory Rate (RR):
  - Titrate to maintain PaCO2 35-45 mmHg, or permissive hypercapnia (PaCO2 50-70 mmHg, pH >7.20) in ARDS to limit VILI.
  - **Rationale:** Higher RR increases VE, reducing PaCO2; lower RR in obstructive diseases (e.g., COPD) prevents dynamic hyperinflation.
- FiO2:
  - Titrate to SpO2 88-95% (PaO2 60-100 mmHg), avoiding hyperoxia (FiO2 >60% for >24h).
  - **Rationale:** Minimizes oxygen toxicity, absorption atelectasis, and ROS-mediated lung injury.

- PEEP:
  - **ARDSNet PEEP-FiO2 Table:** E.g., FiO2 0.5 with PEEP 8-10 cmH2O, FiO2 1.0 with PEEP 14-24 cmH2O.
  - **Optimal PEEP:** Titrate using driving pressure ( $\Delta P = P_{\text{plat}} - \text{PEEP}$ , target <15 cmH2O) or esophageal pressure monitoring (transpulmonary pressure 0-10 cmH2O).
  - **Rationale:** Higher PEEP improves oxygenation by increasing FRC, but risks overdistension ( $P_{\text{plat}} > 30$  cmH2O) and hemodynamic compromise (decreased venous return).
- Inspiratory Time (Ti) and I:E Ratio:
  - Normal: I:E 1:2; ARDS: 1:1 (increase Ti for oxygenation); Obstructive: 1:4-1:5 (prolong Te to prevent air trapping).
  - **Rationale:** Longer Ti increases MAP, improving oxygenation; longer Te reduces auto-PEEP in obstructive diseases.
- Driving Pressure ( $\Delta P$ ):
  - $\Delta P = P_{\text{plat}} - \text{PEEP}$ ; target <15 cmH2O (associated with improved survival in ARDS per Amato et al., NEJM 2015).
  - **Rationale:** Lower  $\Delta P$  minimizes mechanical stress, reducing VILI.
- Auto-PEEP:
  - **Measure with end-expiratory hold:** Auto-PEEP >5 cmH2O indicates air trapping.
  - **Rationale:** Common in obstructive diseases (e.g., COPD, asthma); reduce by lowering RR, prolonging Te, or decreasing VT.

## Troubleshooting Advanced Ventilator Issues

### Elevated Peak Pressure ( $P_{\text{peak}} > 40$ cmH2O):

- Differential:
  - **Airway Resistance:** Bronchospasm, secretions, ETT obstruction ( $R = P_{\text{peak}} - P_{\text{plat}} / \text{Flow}$ ).
  - **Decreased Compliance:** Pneumothorax, ARDS, pulmonary edema ( $C_{\text{st}} = \text{VT} / [P_{\text{plat}} - \text{PEEP}]$ ).
  - **Patient Factors:** Biting ETT, coughing, dyssynchrony.
- Management:
  - **If  $P_{\text{peak}} - P_{\text{plat}} > 10$  cmH2O:** High resistance; suction ETT, administer bronchodilators (albuterol 2.5 mg nebulized), ensure ETT patency.
  - **If  $P_{\text{plat}} > 30$  cmH2O:** Low compliance; reduce VT to 4-6 mL/kg, lower PEEP if overdistension suspected, CXR to rule out pneumothorax.
  - **Dyssynchrony:** Increase sedation (e.g., propofol 20-50  $\mu\text{g/kg/min}$  IV), adjust trigger sensitivity (e.g., flow trigger 1-3 L/min).

- Example: "Patient with Ppeak 48 cmH<sub>2</sub>O, Pplat 34 cmH<sub>2</sub>O (Cst 20 mL/cmH<sub>2</sub>O), CXR negative, VT reduced from 8 to 5 mL/kg, PEEP lowered from 15 to 10 cmH<sub>2</sub>O, Ppeak decreased to 38 cmH<sub>2</sub>O."

## Refractory Hypoxemia (PaO<sub>2</sub>/FiO<sub>2</sub> <150 mmHg Despite FiO<sub>2</sub> >60%) :

---

- Differential:
  - **Inadequate Recruitment:** Low PEEP, atelectasis.
  - **V/Q Mismatch:** ARDS, pneumonia, PE.
  - **Shunt:** Pneumothorax, right-to-left shunt (e.g., PFO).
- Management:
  - **Optimize PEEP:** Use ARDSNet PEEP-FiO<sub>2</sub> table or esophageal pressure-guided PEEP (transpulmonary pressure 0-10 cmH<sub>2</sub>O).
  - **Recruitment Maneuvers:** 40 cmH<sub>2</sub>O for 40 s (if hemodynamically stable).
  - **Prone Positioning:** 16h/day in severe ARDS (PaO<sub>2</sub>/FiO<sub>2</sub> <150, PROSEVA trial, NEJM 2013).
  - **Advanced Modes:** Switch to APRV (P-high 25 cmH<sub>2</sub>O, T-high 5 s) or HFOV (mPaw 30 cmH<sub>2</sub>O, Hz 5).
  - **Treat Underlying Cause:** Antibiotics for pneumonia, thrombolytics for massive PE.
- Example: "Patient with refractory hypoxemia (PaO<sub>2</sub>/FiO<sub>2</sub> 120 on FiO<sub>2</sub> 80%, PEEP 10 cmH<sub>2</sub>O), switched to APRV (P-high 28 cmH<sub>2</sub>O, T-high 5 s, T-low 0.8 s), prone positioning initiated, PaO<sub>2</sub>/FiO<sub>2</sub> improved to 180."

## Hypercapnia with Acidosis (PaCO<sub>2</sub> >50 mmHg, pH <7.35) :

---

- Differential:
  - **Inadequate VE:** Low VT/RR, high dead space (V<sub>d</sub>/VT >0.6).
  - **Dynamic Hyperinflation:** Auto-PEEP in obstructive diseases (COPD, asthma).
  - **Increased CO<sub>2</sub> Production:** Fever, seizures, overfeeding.
- Management:
  - **Increase VE:** Raise RR (e.g., 20 to 24 breaths/min), increase VT if Pplat <30 cmH<sub>2</sub>O.
  - **Reduce Auto-PEEP:** Prolong Te (I:E 1:4-1:5), reduce RR, bronchodilators (albuterol 2.5 mg nebulized q4h).
  - **Decrease CO<sub>2</sub> Production:** Control fever (acetaminophen 650 mg IV), reduce carbohydrate load in TPN.
  - **Permissive Hypercapnia:** Allow PaCO<sub>2</sub> 50-70 mmHg (pH >7.20) in ARDS to minimize VILI.

- Example: "Patient with hypercapnia (PaCO<sub>2</sub> 65 mmHg, pH 7.28), COPD exacerbation, auto-PEEP 8 cmH<sub>2</sub>O, I:E adjusted to 1:4, RR reduced to 12 breaths/min, albuterol given, PaCO<sub>2</sub> decreased to 55 mmHg."

## Patient-Ventilator Dyssynchrony:

---

- Types:
  - **Trigger Dyssynchrony:** Insensitive trigger (e.g., flow trigger >3 L/min), auto-PEEP preventing triggering.
  - **Flow Dyssynchrony:** Inadequate inspiratory flow rate (e.g., <60 L/min in VC).
  - **Cycling Dyssynchrony:** Ti too short/long, leading to double-triggering or breath stacking.
- Management:
  - **Optimize Trigger:** Use flow trigger 1-2 L/min, address auto-PEEP (see above).
  - **Adjust Flow:** Increase flow rate in VC (e.g., 60-80 L/min), or switch to PC mode.
  - **Correct Ti:** Set Ti 0.8-1.2 s, I:E 1:2; in double-triggering, prolong Ti or increase sedation.
  - **Sedation/Analgesia:** Propofol 20-50 µg/kg/min IV, fentanyl 50-100 µg/h IV.
  - **Neuromuscular Blockade:** Cisatracurium 2 µg/kg/min IV (last resort, e.g., severe ARDS with dyssynchrony).
- • Example: "Patient with trigger dyssynchrony (auto-PEEP 10 cmH<sub>2</sub>O, COPD), flow trigger adjusted to 1 L/min, I:E set to 1:4, propofol started, dyssynchrony resolved, auto- PEEP decreased to 4 cmH<sub>2</sub>O."

## Hemodynamic Compromise (Hypotension Post-Ventilation):

---

- Differential:
  - **Decreased Venous Return:** High PEEP/MAP increases intrathoracic pressure.
  - **Tension Pneumothorax:** High PEEP causing barotrauma.
  - **RV Dysfunction:** High PEEP increasing RV afterload (ARDS, PE).
- Management:
  - **Lower PEEP:** Reduce PEEP (e.g., 15 to 10 cmH<sub>2</sub>O) if MAP >20 cmH<sub>2</sub>O.
  - **Fluid Bolus:** NS 500 mL IV to improve preload.
  - **Vasopressors:** Norepinephrine 5-20 µg/min IV if MAP <65 mmHg.
  - **ECHO:** Assess RV function (RV dilation, septal flattening); if present, reduce PEEP, consider inotropes (e.g., milrinone 0.375 µg/kg/min).
  - **CXR:** Rule out tension pneumothorax (chest tube if present).

- Example: "Patient with hypotension (MAP 60 mmHg) post-ventilation (PEEP 15 cmH2O), ECHO showing RV dilation, PEEP reduced to 10 cmH2O, NS 500 mL bolus given, MAP improved to 70 mmHg."

**Table:** Advanced Ventilator Settings by Clinical Scenario

Scenario	Mode	VT (mL/kg IBW)	RR (breaths/min)	FiO2 (%)	PEEP (cmH2O)	I:E Ratio	Rationale
Severe ARDS (PaO2/FiO2 <150)	AC-PC or APRV	4-6	12-18	60-100	1:1 (AC), T-high 5-6 s (APRV)	Low VT, high PEEP to minimize VILI, maximize recruitment; APRV for refractory hypoxemia	
AECOPD with Hypercapnia	AC-VC	6-8	8-12	40-60	5-8	1:4-1:5	Low RR, prolonged Te to reduce auto-PEEP, avoid high PEEP to prevent hyperinflation
Status Asthmaticus	AC-VC	6-8	6-10	40-60	0-5	1:4-1:5	Very low RR, minimal PEEP to avoid air trapping, maximize exhalation
Neuromuscular Failure (e.g., GBS)	AC-VC	6-8	12-16	30-50	5-8	1:2	Full support, normal I:E, moderate RR to prevent hypercapnia
Weaning (SBT)	PSV	N/A	N/A	30-40	5	N/A	PS 5-15 cmH2O to support spontaneous breaths, RSBI <105 predicts success



Scenario	Mode	VT (mL/ kg IBW)	RR (breaths/ min)	FiO2 (%)	PEEP (cmH2O)	I:E Ratio	Rationale
Refractory Hypoxemia (Rescue)	HFOV	1-3	5-8 Hz	60-100	mPaw 30-35	N/A	Ultra-low VT, high frequency to minimize VILI, improve oxygenation

## Key Pearls for Advanced Practice

### • Physiology:

- Oxygenation: Increasing O2
  - Optimize FiO2, PEEP, MAP; target SpO2 88-95% to avoid toxicity.
- Ventilation: Decreasing CO2
  - Adjust VE (VT × RR) to maintain PaCO2 35-45 mmHg; use permissive hypercapnia in ARDS (pH >7.20).
- Lung Mechanics:
  - Monitor Pplat (<30 cmH2O), ΔP (<15 cmH2O), auto-PEEP to prevent VILI.

### • Settings:

- ARDS:
  - Low VT (4-6 mL/kg), high PEEP (12-18 cmH2O), PC or APRV for pressure control.
- Obstructive ex COPD:
  - Low RR (6-12 breaths/min), I:E 1:4-1:5, minimal PEEP to reduce air trapping.
- Weaning:
  - PSV with RSBI <105, FiO2 30-40%, PEEP 5 cmH2O.

### • Troubleshooting:

- High Ppeak:
  - Differentiate resistance (Ppeak - Pplat) vs. compliance (Pplat); address obstruction, reduce VT/PEEP.
- Hypoxemia:
  - Optimize PEEP (esophageal pressure-guided), prone positioning, APRV/HFOV for rescue.
- Hypercapnia:
  - Increase VE, reduce auto-PEEP with prolonged Te, treat CO2 production.

- Dyssynchrony:
  - Adjust trigger, flow,  $T_i$ ; use sedation or paralysis as needed.

## • **Monitoring:**

---

- ABG q4-6h: Assess  $PaO_2$ ,  $PaCO_2$ , pH.
- Ventilator Parameters:  $P_{peak}$ ,  $P_{plat}$ , auto-PEEP,  $\Delta P$ .
- Hemodynamics: ECHO for RV strain, MAP >65 mmHg.

## Advanced Clinical Scenarios

### Scenario 1: Severe ARDS with Refractory Hypoxemia

---

- Presentation: A 55-year-old male with ARDS secondary to sepsis ( $PaO_2/FiO_2$  120 on  $FiO_2$  80%, PEEP 12 cmH<sub>2</sub>O) remains hypoxemic despite conventional ventilation.
- Ventilator Strategy:
  - **Switch to APRV:** P-high 28 cmH<sub>2</sub>O, P-low 0 cmH<sub>2</sub>O, T-high 5 s, T-low 0.8 s (set to 75% peak expiratory flow).
  - $FiO_2$  90%, titrate to  $SpO_2$  88-92%.
  - Prone positioning initiated for 16h/day. Driving pressure monitored ( $\Delta P$  14 cmH<sub>2</sub>O).
- Rationale: APRV maximizes recruitment, minimizes cyclic collapse; proning improves V/Q matching (PROSEVA trial).
- Outcome:  $PaO_2/FiO_2$  improves to 180 after 12h,  $FiO_2$  reduced to 60%, patient stabilized.

### Scenario 2: Status Asthmaticus with Dynamic Hyperinflation

---

- Presentation: A 30-year-old female with status asthmaticus ( $PaCO_2$  70 mmHg, pH 7.25, auto-PEEP 12 cmH<sub>2</sub>O) on AC-VC (VT 8 mL/kg, RR 16 breaths/min, PEEP 5 cmH<sub>2</sub>O).
- Ventilator Strategy:
  - VT reduced to 6 mL/kg (420 mL for 70 kg IBW).
  - RR decreased to 8 breaths/min, I:E set to 1:5 ( $T_i$  0.5 s,  $T_e$  2.5 s).
  - PEEP maintained at 5 cmH<sub>2</sub>O.
  - Albuterol 2.5 mg nebulized q2h, IV steroids (methylprednisolone 125 mg IV).
- Rationale: Low RR and prolonged  $T_e$  reduce air trapping, auto-PEEP; bronchodilators address airway resistance.
- Outcome: Auto-PEEP decreases to 6 cmH<sub>2</sub>O,  $PaCO_2$  to 50 mmHg, pH 7.35 after 4h, patient stabilizes.

### Scenario 3: Weaning Failure in Neuromuscular Disease

---

- Presentation: A 45-year-old male with Guillain-Barré syndrome on AC-VC (VT 7 mL/kg, RR 14 breaths/min, FiO<sub>2</sub> 40%, PEEP 5 cmH<sub>2</sub>O) fails SBT (RSBI 120, RR 30 breaths/min, VT 200 mL).
- Ventilator Strategy:
  - **Switch to SIMV:** VT 7 mL/kg, RR 10 breaths/min, PS 10 cmH<sub>2</sub>O, FiO<sub>2</sub> 40%, PEEP 5 cmH<sub>2</sub>O.
  - Increase sedation (fentanyl 50 µg/h IV) to reduce RR.
  - Monitor RSBI daily, aim for <105.
- Rationale: SIMV provides partial support, reduces WOB; PS augments spontaneous breaths, preventing fatigue.
- Outcome: After 48h, RSBI decreases to 90, patient passes SBT (RR 18 breaths/min, VT 350 mL), extubated successfully.

**Visit: [webcheatsheets.com](https://www.webcheatsheets.com) for more education, fun resources and 10 category 1 AAPA CME credit!**

© Hospital Medicine Cheat Sheets ([medcheatsheets.com](https://www.medcheatsheets.com)). For educational purposes only. Do not redistribute or sell. Neither the author nor the company is liable for real-world implications. AI was used in development